



Heartland
Medical Sales & Services, LLC

www.heartlandmedical.com

Heartland LS 10 LED 10-Watt LED Light Source



The Heartland LS 10 LED 10-Watt Light Source provides beautiful, long-lasting white light to your headlights, retractors and lighted instruments. Lightweight, small and unobtrusive, yet powerful. The LS 10 LED can sit table-top or can be attached to any steel surface with its base-mounted earth magnet.

Lamp: To 4,350 Foot Candles*
Color Temp: 6,000° Kelvin
Lamp Life: 25,000 Hours
Dimensions: 3.25" x 2.5" x 1.75" (82mm x 64mm x 44mm)
Weight: Less than 10 oz. / .28kg
Power: 100v - 240v 50/60 Hz