



Heartland
Medical Sales & Services, LLC

www.heartlandmedical.com

Heartland PS 4500 LED 3-Watt Portable LED Headlight System



The Heartland PS 4500 LED Portable LED headlight system brings brilliant, even, adjustable LED intensity into the field. With variable intensity control, Pinpoint Spot Adjustment and a sterilizable positioning bar, the PS 4500 LED headlight brings versatility and functionality to the operating room.

Lightweight, comfortable headband options are available as well as a transformer option.

Lamp: 3-Watt LED / 4,500 Foot Candles
Color Temp: 6,000° Kelvin (3,200° Kelvin optional)
Lamp Life: 50,000 Hours
Battery Life: 7.5 to 100 hours
Power: 100v - 240v Charger